



Annual Evacu B Maintenance Checklist

Date:

Institution Name / Address:

Check performed by:

Chair Serial #:

Check if in working condition:

- Track drops and locks
- White belt turn when gliding down the stairs
- Interior pocket walls are attached
- Pocket tethers are loose
- Wrist strap is intact
- Yellow decals
- Yellow instruction card on brake handle
- No rips in the red apron material
- All Velcro for headrests are firmly attached – to headrest and apron
- All headrests are firmly attached
- Brake locks into place when handle is released
- Casters are secured tightly

January 20, 2017

Signature: _____ Date: _____