



Wall Bracket Instructions

- Mount the EVACU B to the wall using the provided wall bracket
- Measuring from the bottom of the bracket, mount the wall bracket at least 62" up from the floor
- Make sure the wall bracket is mounted at least 5" away from a corner wall or another wall bracket
- When hung, the EVACU B occupies a space that is 68" Tall, 24" Wide and will come out 12" off of the wall



EVACU B Specifications:

Weight: 26 lbs
Length: 59"
Width: 23 ¾"
Height: 12"

*Dimensions from a collapsed position

Doug Gervais

President/Doro Inc.

August 22, 2016